

## ***Reactivity to novel Influenza A(H1N1) virus***

2009.10

*Following 30 strains of novel Influenza A (H1N1) virus has been confirmed to be positive by this test kit.*

<b>Subtype</b>	<b>Name of strain</b>	<b>Subtype</b>	<b>Name of strain</b>
Pandemic H1N1	A/Osaka/55/2009(H1N1)pdm	Pandemic H1N1	A/Osaka/125/2009(H1N1)pdm
Pandemic H1N1	A/Osaka/68/2009(H1N1)pdm	Pandemic H1N1	A/Osaka/126/2009(H1N1)pdm
Pandemic H1N1	A/Osaka/71/2009(H1N1)pdm	Pandemic H1N1	A/Osaka/143/2009(H1N1)pdm
Pandemic H1N1	A/Osaka/72/2009(H1N1)pdm	Pandemic H1N1	A/Osaka/144/2009(H1N1)pdm
Pandemic H1N1	A/Osaka/78/2009(H1N1)pdm	Pandemic H1N1	A/Osaka/146/2009(H1N1)pdm
Pandemic H1N1	A/Osaka/84/2009(H1N1)pdm	Pandemic H1N1	A/Osaka/148/2009(H1N1)pdm
Pandemic H1N1	A/Osaka/101/2009(H1N1)pdm	Pandemic H1N1	A/Osaka/157/2009(H1N1)pdm
Pandemic H1N1	A/Osaka/107/2009(H1N1)pdm	Pandemic H1N1	A/Osaka/164/2009(H1N1)pdm
Pandemic H1N1	A/Osaka/110/2009(H1N1)pdm	Pandemic H1N1	A/Osaka/167/2009(H1N1)pdm
Pandemic H1N1	A/Osaka/112/2009(H1N1)pdm	Pandemic H1N1	A/Osaka/168/2009(H1N1)pdm
Pandemic H1N1	A/Osaka/114/2009(H1N1)pdm	Pandemic H1N1	A/Osaka/171/2009(H1N1)pdm
Pandemic H1N1	A/Osaka/115/2009(H1N1)pdm	Pandemic H1N1	A/Osaka/172/2009(H1N1)pdm
Pandemic H1N1	A/Osaka/116/2009(H1N1)pdm	Pandemic H1N1	A/Osaka/174/2009(H1N1)pdm
Pandemic H1N1	A/Osaka/118/2009(H1N1)pdm	Pandemic H1N1	A/Osaka/176/2009(H1N1)pdm
Pandemic H1N1	A/Osaka/119/2009(H1N1)pdm	Pandemic H1N1	A/Osaka/193/2009(H1N1)pdm

**Tested at OSAKA PREFECTURAL INSTITUTE OF PUBLIC HEALTH**

**MIZUHO MEDY Co.,Ltd.**